

# LAB PRODUCT CATALOGUE



*Global Standards in Dental Technology,  
Trusted by Professionals .*



11+  
Offices



31+ Global  
Brands



200+  
Products

EXCLUSIVE  
DISTRIBUTOR

# Profeta<sup>®</sup> VP 100 PRO

SLM Technology 3D Metal Printer

## VERSATILE PRINTER 100 PRO



# Dual mode machine for easy powder switching

## VERSATILE PRINTER 100 PRO (VP100 PRO)

Building Volume	Φ120mm, z=70mm
Layer Thickness	20-50 μm
Scanning Speed	15m/s(max)
Productivity	130copings/4-5 hours
	10RPD frames/3.5-4 hours
Laser Type	Yb-fiber laser, 200W/300W
Laser Beam Diameter	40~60 μm
Power	AC220V, 12A
	2.5KW
Shielding Gas	Argon/Nitrogen
Machine Dimension	900*680*1760mm
Weight	480kg
Operating Temperature	18~30 °C
Materials	CoCr, Ti, TC4
Applicable Products	Copings, RPD frames, implant bar surgical guide, orthodontic appliance etc.

### DEDICATED DENTAL METAL 3D PRINTER

- ◆ **Dual-mode** printer supports Cobalt-Chromium and Titanium, and the switching time is only 45 minutes;
- ◆ Easy to start with Cobalt-Chromium mode and be able to upgrade to dual-mode later on;
- ◆ Equipped with glove box and built-in cleaner to implement **closed powder cleaning** and to avoid harmful dust;
- ◆ High printing quality for RPD frames, copings, post & core, surgical guide, orthodontic appliance etc;
- ◆ Proved stability in domestic and global market.



### ONE-STOP SOLUTION PACKAGE

- ◆ The one-stop solution package includes design software, slicing software, furnace, automatic screening cabinet and nitrogen generator;
- ◆ Our solution achieves synergy effects for better quality and efficiency.



### CUSTOMER VALUE

- ◆ Designed for **small-middle scale dental labs**, and greatly supports all scale labs which require frequent powder switching;
- ◆ Compact machine, requires minimum facility area, and easy to transport through elevator;
- ◆ One dual-mode printer serves the purpose of two machines;
- ◆ Closed powder cleaning function is designed for better HSE.



EXCLUSIVE  
DISTRIBUTOR

 Profeta<sup>®</sup> IE 150

SLM Technology 3D Metal Printer

INTELLIGENT AND EXPERT  
EFFICIENT

SINGLE LASER



Also Available in  
DUAL LASER

# Single lasers for CoCr or Ti

## INTELLIGENT & EXPERT PRINTER (IE 150)

Building Volume	Φ158 mm, z=70 mm
Layer Thickness	20-50 μm
Scanning Speed	15m/s(max)
Productivity	220 copings/3.5-4.5 hours
	15 RPD frames/3-3.5 hours
Laser Type	Yb-fiber laser*2, 200W/300W
Laser Beam Diameter	40~60 μm
Power	AC220V, 16A
	3.5 KW
Shielding Gas	Argon/Nitrogen
Machine Dimension	1070*670*1750 mm
Weight	500 kg
Operating Temperature	18~30 °C
Materials	CoCr, Ti, TC4
Applicable Products	Copings, RPD frames, implant bar, surgical guide, orthodontic appliance etc.

### DEDICATED DENTAL METAL 3D PRINTER

- ◆ Equipped with dual lasers for faster printing speed
- ◆ Larger building volume (Φ158 mm, z=70 mm) for higher production
- ◆ Higher powder utilization
- ◆ Support two printing layer thickness (20-50μm)
- ◆ Proved stability in domestic and global market.



### ONE-STOP SOLUTION PACKAGE

- ◆ The one-stop solution package includes design software, slicing software, furnace and nitrogen generator;
- ◆ Our solution achieves synergy effects for better quality and efficiency.



### CUSTOMER VALUE

- ◆ Designed for all scale dental labs
- ◆ Higher powder utilization, lower production cost per unit;
- ◆ Proved machine stability.



EXCLUSIVE  
DISTRIBUTOR

 Profeta<sup>®</sup> C 150

Dedicated Dental Metal 3D Printer

**CLASSICAL DENTAL PRINTER**



# Single Lasers for CoCr

## CLASSICAL DENTAL PRINTER (C150)

Building Volume	Φ158 mm, z=70 mm
Layer Thickness	20-50 μm
Scanning Speed	15m/s(max)
Productivity	220 copings/4-4.5 hours
	15 RPD frames/3.5-4 hours
Laser Type	Yb-fiber laser (Single/Dual) 300W/500W
Laser Beam Diameter	40~60 μm
Power	AC220V, 16A
	3.5 KW
Shielding Gas	Nitrogen
Machine Dimension	1070*670*1750 mm
Weight	500 kg
Operating Temperature	18~30 °C
Materials	CoCr
Applicable Products	Copings, RPD frames, implant bar, surgical guide, orthodontic appliance etc.

### DEDICATED DENTAL METAL 3D PRINTER

- ◆ Equipped with dual lasers for faster printing speed;
- ◆ Larger building volume (Φ158 mm, z=70 mm) for higher production
- ◆ Higher powder utilization
- ◆ Support two printing layer thickness (20-50μm);
- ◆ Proved stability in domestic and global market.



### ONE-STOP SOLUTION PACKAGE

- ◆ The one-stop solution package includes design software, slicing software, furnace and nitrogen generator;
- ◆ Our solution achieves synergy effects for better quality and efficiency.



### CUSTOMER VALUE

- ◆ Designed for all scale dental labs
- ◆ Higher powder utilization, lower production cost per unit;
- ◆ Proved machine stability.



# Profeta<sup>®</sup> C 120

## High-Performance Compact Model



## PREMIUM DENTAL PRINTER

### PROFETA C120

Building Volume	Φ120mm , z=70mm
Layer Thickness	20-50 μm
Scanning Speed	15m/s(max)
Laser Type	Yb-fiber laser, 300W
Laser Beam Diameter	40-60 μm
Power	AC 220V, 12A
	2.5 KW
Shielding Gas	Nitrogen
Machine Dimension	900*680*1760 mm
Weight	380 kg
Operating Temperature	18-30°C
Materials	CoCr

# AUXILIARY EQUIPMENT



Name	Sieving Machine
Purpose	Used for powder recycling and filtering out large metal particles generated during the printing process.
Power	220 V, 50 Hz, 50w
Weight	10kg
Machine Dimension	320*330*380mm



Name	Box-type Resistance Furnace
Purpose	Used for stress relief of 3D printed components.
Power	220V, 5kw
Weight	90kg
Internal Dimension	300*200*120mm
Max. Temperature	1200°C



Name	Bandsaw
Purpose	Cut the printed part from the build plate.
Power	220V, 50Hz, 680w
Weight	15.5kg
Machine Dimension	610*340*410mm
Saw Band Size	1332*13*0.65mm




Name	Nitrogen Generator (Optional)
Purpose	To isolate oxygen during the printing process, prevent oxidation during laser melting, and ensure printing quality.
Power	220V, 50Hz, 50w
Weight	20kg
Machine Dimension	550*250*800mm
Nitrogen Output	10 L/min (adjustable)
Nitrogen Purity	99.99% (adjustable)

# WESHINING

DENTAL AUTOMATED POLISHING SYSTEM

## Plasma Polishing Machine EPC-30



- Compact
- Cost-effective
- High-quality
- Time Saving

**Single Mode Machine  
either CoCr or Ti**

### PRODUCT SPECIFICATIONS

Model	EPC-30
Power	AC220V
Capacity	Crown: 3-18, RPD: 1-2
Machine Size	400X600X930 mm
Production	Small-Medium Production
Industry	Dentistry, Medical

# One Step Solution:




Surface Treatment Services



Machine Development




Machine Testing



After-sales Training Service

## Fast Speed:

Reduced reprocessing time and costs



### DENTAL POLISHING SYSTEM:

Takes only 7-10 mins




Before Polishing



After Polishing

### DIGITAL DENTAL CASES:



Before



After

Powder (800g):	RPD	Replacement of Powder
CoCr	250 RPD	1000-1200 Crown
Ti	80 RPD	300-400 Crown

# INTRAORAL SCANNER 3DS V5



C-Type Data Cable



Extended Scanning Tip



Metal Scan



Quicker Scanning Speed



Case Cloud Sharing




Ultra Lightweight

## PRODUCT SPECIFICATIONS

Handpiece Size	243x40x32mm
Connection Method	Wired Type-C connection
Handpiece Weight	160g
Arcbolames Per Second	20 FPS
Scanning Accuracy	Single Crown $\leq$ 10 $\mu$ m, Full Arch $\leq$ 20 $\mu$ m
Scanner Tips	3 standard /1 mini /1 Cover
Scanning Area	Standard 16x14mm/Mini 12x10mm
Imaging Method	Optical Video Capture
Quick Reset	Yes
Model Optimization	Yes
Anti-fog Method	Fan-cooling
Calibration Tip	1 Set (Optional)

# INTRAORAL SCANNER 3DS V3 PRO

## Select your case workflow



Calibration Free Scanner



Fashion & Light Handpiece



Scan Body Smart Matching



AI Fast Scan



Humanized Scanner Tip



Metal Scan

**500+** Installations within a year

## PRODUCT SPECIFICATIONS

Handpiece Size	240x49.5x30.5mm
Connection Method	Direction Connect Plug & Play
Handpiece Weight	210g
Scanning Depth	15 mm
Scanner Light Source	Advanced LED
Small Tips Size	80 X 16 X 11 mm
Operating Time	>30000 hrs
Imaging Method	Synchronous 3D Video
Remote Control Function	Yes
Data Optimization	Yes
Export File	STL, PLY, OBJ
Calibration	Not Required

# 4-Axis Wet Dental Mill

## DWX-43W




Mill Up to 6  
Pin-Type Blocks



Efficient Operation



DGSHAPE CAM  
Software Included

**Ti**

Mill Titanium (only pre-  
milled abutments)



Easy to Use



Remote Monitor

## PRODUCT SPECIFICATIONS

Millable material *1	Glass ceramics, Titanium alloy, Composite resins, PMMA, PEEK, Sintered zirconia
Loadable workpiece shape	Pin type (Width x depth x height): Max. 40 x 20 x 20 mm (1.5 x 0.8 x 0.8 in.)
Operating speed	XYZ: 6 to 1,800 mm/min (0.24 to 70.9 in./min.)
Attachable milling bur	Shank diameter: 3 mm, dedicated milling bur
Compatible compressed air	0.18 to 0.22 Mpa
Spindle rotating speed	15,000 to 60,000 rpm
Rotary axis travel angle	A:±36°
Number of milling burs housed	6

# 5-Axis Dry Dental Mill

## DWX-53D



5-Axis Dental Milling



Up To 15 Burs Auto  
Tool Changer



Built-in Temperature  
Sensor



Open CAD/CAM  
System



User-Replacement  
Spindle



Remote Monitor

## PRODUCT SPECIFICATIONS

Millable Material		Zirconia, Wax, PMMA, Composite resin, PEEK, PEKK Gypsum, Glass fibre re-inforced resin, CoCr sinter metal, Polycarbonate
Loadable workpiece Shape	Disc-type (with step)	Step section: 98.5 mm (3.88 in.) (diameter), $10 \pm 0.3$ mm ( $0.39 \pm 0.01$ in.) (height) Body section: 95 mm (3.74 in.) (diameter), 60 mm (2.36 in.) max. (height)
	Pin-type	Pin diameter: 6 mm (0.24 in.)
Operating speed		X, Y: 3,000 mm/min (118.11 in./min) max., Z: 2,700 mm/min (106.29 in./min) max.
Spindle motor		Brushless DC Motor
Spindle rotating speed		6,000 to 30,000
Rotary axis travel angle		A: $\pm 360^\circ$ , B: $\pm 35^\circ$
No. of milling burs housed		15

# 5-Axis Dry Dental Mill

## DWX-52D PLUS



5-Axis Dental Milling



Up To 15 Burs Auto  
Tool Changer



**90°**

90° Milling (with the  
Open Edge Adapter:  
ZMA2-52D)



Open CAD/CAM  
System



User-Replacement  
Spindle



Remote Monitor

### PRODUCT SPECIFICATIONS

Millable Material		Zirconia, Wax, PMMA, Composite resin, PEEK, PEKK, Gypsum, Glass fiber reinforced resin, CoCr sinter metal, Polycarbonate
Loadable workpiece Shape	Disc-type (with step)	Step section: 98.5mm (3.88 in.) (diameter), Height: 10mm (0.39 in.), Body section: 95 mm (3.74 in.) (diameter), Height max. 60 mm (2.36 in.)
	Pin-type	Pin diameter: 6 mm (0.24 in.)
Operating speed		XYZ axis: 6 to 1800 mm/min. (0.24 to 70.9 in/min.)
Spindle motor		Brushless DC Motor
Spindle rotating speed		6,000 to 30,000
Rotary axis travel angle		A:±360°, B: ±30°
No. of milling burs housed		15

# 5-Axis Dental Mill with Disc Changer

## DWX-52DI PLUS



5-Axis Dental Milling



Up To 8 Pin-Type  
Blocks



**90°**  
90° Milling  
Range



Open CAD/CAM  
System



C-type Clamping  
System



Entry-Level Price

### PRODUCT SPECIFICATIONS

Millable material *1		Zirconia, Wax, PMMA, Composite resin, PEEK, PEKK Gypsum, Glass fibre re-inforced resin, CoCr sinter metal, Polycarbonate
Loadable workpiece Shape	Disc-type (with step)	Step section: 98.5 mm (3.88 in.) (diameter), height 10 mm (0.39 in.), Body section: 95 mm (3.47 in.) (diameter), height max. 60 mm (2.36 in.)
	Pin-type	Pin diameter: 6 mm (0.24 in.) *The optional pin-type material adapter ZV-52D is required.
Operating speed		XYZ axis: 6 to 1800 mm/min. (0.24 to 70.9 in./min.)
Spindle motor		Brushless DC Motor
Spindle rotating speed		6,000 to 30,000
Rotary axis travel angle		A: ±360°, B: ±30°
No. of milling burs housed		15

# SHINING 3D DENTAL

## Autoscan-DS-EX PRO (H)



Unparalleled  
Details



Ultra-High  
Accuracy



Outstanding  
Scan Speed



Flexible  
Workflow



### TECHNICAL SPECIFICATIONS

Camera Resolution	2.5 MP
Accuracy	8µm
Scan Time	Bite scan: 6.5s; Upper/Lower jaw scan: 14s; 1-4 dies scan: 21.5s; 5-8 dies scan: 25.5s; Impression scan: 31s;
Articulator Scan	Support articulator transfer
Output Format	STL, PLY, OBJ
Light Source	Blue Light
Dimensions	420*260*270 mm
Weight	5.0 Kg
Interface	USB 3.0
Power Supply	DC 24V

# SHINING 3D DENTAL

## Autoscan-DS-EX PRO (C)

HIGH EFFICIENCY, COLOR TEXTURE EASY FOR RPD DESIGN



Portable



Outstanding  
Scan Speed



Compatible



### 2 STEP SCAN



1. Occlusion



2. All In one

### TECHNICAL SPECIFICATIONS

Camera Resolution	2.4 MP
Accuracy	< 10µm
Scan Time	5 s for bite scan; 9.5 s for upper / lower jaw scan; 15 s for 1-4 dies scan; 17 s for 5-8 dies scan; 22 s for impression scan.
Output Format	STL, OBJ, PLY
Light Source	Blue Light
Dimensions	260*270*420mm
Weight	5.0 Kg
Interface	USB 3.0
Power Supply	DC 24V

# exocad

## DentalCAD 3.2 Elefsina

EXPERIENCE THE NEXT LEVEL OF DIGITAL DENTISTRY



For AI services upgrade to the Feature Release now!



**ADVANCED AUTOMATION**



**ELEVATED SMILE CREATOR**



**PRECISION TOOLS**



**SEAMLESS WORKFLOWS**

### TECHNICAL SPECIFICATIONS

Modules	Smile Creator, Bite Splint, Full Denture, Post & Core
Output Format	STL, OBJ, PLY
Photo Matching	Advanced 3D-photo alignment
Add-On Compatibility	Ivotion™, Milling Integration, iTero, exoplan
Integration	Compatible with all major intraoral scanners
Automation	Auto Articulator, Scan Margin Repair
System Compatibility	Works with all major CAD/CAM dental systems

## 2 Packages Available for DentalCAD

Product		Core Lab Version	Advanced Lab Version	Implant Lab Version	Ultimate Lab Version
	DentalCAD Core Version	✓	✓	✓	✓
	Virtual Articulator	Add-on	✓	✓	✓
	Provisional Module	Add-on	✓	✓	✓
	TruSmile	Add-on	✓	✓	✓
	ZRS Tooth Library	Add-on	✓	✓	✓
	Implant Module	Add-on	Add-on	✓	✓
	Bar Module <sup>2</sup>	Add-on	Add-on	✓	✓
	DICOM Viewer	Add-on	Add-on	✓	✓
	Model Creator <sup>4</sup>	Add-on	Add-on	Add-on	✓
	Smile Creators <sup>5</sup>	Add-on	Add-on	Add-on	✓
	FullDenture Module	Add-on	Add-on	Add-on	✓
	PartialCAD Module	Add-on	Add-on	Add-on	✓
	Bite Splint Module	Add-on	Add-on	Add-on	✓
	Jaw Motion Import <sup>6</sup>	Add-on	Add-on	Add-on	✓
	In-CAD Nesting	Add-on	Add-on	Add-on	Add-on
	Nesting (includes in-CAD Nesting)	Add-on	Add-on	Add-on	Add-on



## X-Mill 500 Plus Milling Machine

A newly developed 5-axis automatic engraving and milling equipment. The high-sealing processing chamber has good dust-proof performance and wide applicability. It can process more materials and quickly cut hard materials such as PMMA/PEEK; it has stronger situational awareness and intelligent monitoring; the door opens and stops automatically, making it safe and worry-free. It is not only a change in appearance design, but also more intelligent, safer, wider applicability, and significantly improved performance.

- ◆ Outstanding Performance
- ◆ Intelligent monitoring
- ◆ Safe operating
- ◆ Speed processing for hard materials
- ◆ New appearance design



### TECHNICAL SPECIFICATIONS

Name	X-mill 500 plus
Dimension	(L*W*H): 630*730*700(mm)
Net weight	Approximately 197 kg
Milling Range	XYZ:110/160/80mm A:±30° B:360°
Max. Speed	60,000 rpm
Input Voltage	Single Phase AC 220V
Axis Quantity	5
Milling Mode	Dry milling
Max feeding Speed	XYZ: 6000 mm/min
Bur Changing	Automatic
Bur capacity	8
Spindle Cooling	Water cooling
Materials	zirconia, wax, acril, alluminium , soft metal
Spindle Power	1.8 kw

## 500 SE Milling Machine

Dimension	585X530X695(mm)
Net weight	90 kg
Milling Range	XYZ:170/120/85mm A:±30° B:360°
Max. Speed	60,000 rpm
Input Voltage	AC 220V
No. of Axis	5
Milling Mode	Dry milling
Max Speed	XYZ: 3000 mm/min
Bur Changing	Automatic
Bur capacity	6
Spindle Cooling	Air cooling
Materials	Zirconia, Wax, Aluminium, Soft metal
Spindle Power	0.5 kw



## D5 Milling Machine

Dimension	624*710*720(mm)
Net weight	150 kg
Milling Range	XYZ: 130/100/94mm, A: 360° B: ±30°/-95°
Max. Speed	60,000 rpm
Input Voltage	AC 220V
No. of Axis	5
Milling Mode	Wet and dry
Max Speed	XYZ: 6000 mm/min
Bur Changing	Automatic
Bur capacity	12
Spindle Cooling	Air cooling
Materials	Zirconia, Premill-Abutment, Glass Ceramic
Spindle Power	0.5 kw



## 600 Milling Center

Dimension	1110x820x1820(mm)
Net weight	850kg
Milling Range	XYZ: 320/130/150mm A: 360° B: ±36°
Max. Speed	60,000 rpm
Input Voltage	AC 220V
No. of Axis	5
Milling Mode	Wet milling
Max Feeding Speed	XYZ: 6000mm/min
Bur Changing	Automatic
Bur capacity	12
Spindle Cooling	Water cooling
Materials	CoCr/NEM, Titanium, Aluminium, Soft Metal
Spindle Power	2.5kw



## 630 Milling Centre

Dimension	800X650X1750(mm)
Net weight	450kg
Milling Range	XYZ:240/180/140mm A : 360° B : ±30°
Max. Speed	60,000 rpm
Input Voltage	AC 220V
No. of Axis	5
Milling Mode	Wet milling
Max Speed	XYZ: 6000 mm/min
Bur Changing	Automatic
Bur capacity	24
Spindle Cooling	Water cooling
Materials	Titanium, PMMA, Wax, Glass Ceramic, Premill- Abutment
Milling time	Glass ceramic: 17mins, Wax: 3mins, Premill abutment: 25mins
Spindle Power	2.5 kw



# CHAMPION Equipments

## CHAMPION SINTERING FURNACES

### MODEL: ZTCF-30B PLUS-FAST SINTERING



Fast Sinter	: 65 mins (including heating, keeping warm, and cooling)
Conventional sinter	: 3 layers of crucible, 150 units
Max. Temperature	: 1730°C
Temperature control accuracy	: ±1°C
Touch Screen	: 7" TFT
Max power consumption	: 3 KW
Power Supply	: 220 V 50/60 HZ & 110 V, 50/60 HZ
Heating element	: High purity heater specially for dentistry
Capacity	: Ø 110 ×120 mm
Conventional sintering program	: No.1
Rapid sintering program	: No.2
Cleaning program	: No.40
Program	: 40 programs can be stored, with power-off memory function
Size (L*W*H)	: 590*470*850 mm (L*W*H)
Net weight	: 85 KG

## FAST SINTERING FURNACE

Dental Lab Equipment Can Sinter Zirconia Glass Ceramics And Glaze, it is a high-end sintering furnace, which can be used for sintering zirconia ceramics, glass ceramics and glazes.

### MODEL: F2 Fast Sintering Furnance

Max. Temperature	: 1650°C
Max. Heating Rate	: 100°C/Min
Heating Element	: Silicon Molybdenum Bar
Temperature Control accuracy	: ±1°C
Temperature Controller	: Separate touch screen, 20-segment programmable, Very easy and simple operation
Operation	: Manual
Furnace's size	: 420*320*600mm (smaller & lighter)
Rated voltage	: Interflow 220 VAC±20%
Rated Power	: 2500w



## CHAMPION PORCELAIN FURNACES

### MODEL: ZTCF-20B, CERAMIC FURNACE WITH VACUUM PUMP



Porcelain Furnace temperature control is accurate, the internal temperature of the furnace is uniform, the automatic lifting tray design has good sealing performance, and it is equipped with a low-noise oil-less vacuum pump. Porcelain Furnace is suitable for all kinds of ceramic powder on the market, such as SHOFU, IVOCLAR, VITA, NORITAKE, etc.

Max. Temperature	: 1300°C
Max. Vacuum	: -98 Kpa
Heating wire	: Long lifetime heating wire
Heating rate	: ≤140°C/min
Temperature control accuracy	: ±1°C



## CHAIRSIDE PORCELAIN FURNACE

### ZTCF-20B PRO

Max. Temperature: 1300°C  
 Max. Vacuum: -98Kpa  
 Heating Rate: ≤140°C/min  
 Temperature Accuracy: ±1°C  
 Power: 1600W  
 Voltage: 220V, 50/60HZ  
 Chamber Size: Ø100mm\*70mm  
 Size: 460\*220\*300 mm ( L\*W\*H )  
 Total Program No.: 100 prog. (extend to 500 prog.)  
 Net Weight: 18 KG



## CHAMPION SUCTION



Frequency	50Hz
Size	43.5*41*60cm
Specification	50*50*80cm
Drive	Motor
Voltage	220V
Power	130W
Weight	29kg
Production Capacity	2000units/Year
Certification	CE, CB, ETL, RoHS, UR

### Details

- ❑ 130W Heavy-Duty Motor
- ❑ High Efficiency Filter Element
- ❑ High Quality Sponge
- ❑ Super Strong Wind Speed.

## OIL FREE/LOW NOISE AIR COMPRESSOR

### Model: MEDI+Air-850

Noise : 55 db  
 Tank : 40 L  
 Capacity : 156 L/Min  
 V/Hz : 220V/50Hz  
 Power : 850 W  
 Weight : 31 kg  
 Size : 440 x 440 x 570 L/Min  
 Pressure : 6-8 bar  
 HP : 1.10



## COMPRESSOR WITH DESICAAANT AIRDRYER

### Model: RP03



- ❑ Fitted with 12 cfm desiccant Air dryer and all standard Accessories.
- ❑ Available in 60 Liters tank size.
- ❑ Flow rate 300LPM Maximum
- ❑ Pressure 8 Bar
- ❑ Power 3HP 230V-50Hz (simultaneous on)

## Saturn 4 Ultra (16K)

### Specification

- ❑ XY Resolution: 14×19 μm (15120 × 6230)
- ❑ Build Volume: 211.68 × 118.37 × 220 mm
- ❑ Smart Tank Heating at 30 °C
- ❑ High Speed with Tilt Release Technology
- ❑ Auto-Leveling, Plug-N-Play
- ❑ 4.0-inch Capacitive Touch Screen
- ❑ Upgraded AI Camera with Light for better monitoring and Timelapse Recording
- ❑ MAX 150mm/h Printing Speed



## Saturn 4 (12K)

### Specification

- ❑ XY Resolution: 11520\*5120 (19×24 μm)
- ❑ Build Volume: 218.88×122.88×260 mm
- ❑ COB + Fresnel Collimating Lens
- ❑ Auto Leveling & Intelligent Mechanical Sensor
- ❑ User-Friendly Structure & Advanced Cooling System
- ❑ Large Print Size of 218.88×122.88×220 mm<sup>3</sup>
- ❑ 3.5-inch Resistive Touch Screen
- ❑ 9H hardness tempered glass
- ❑ MAX 70mm/h Printing Speed



## Mercury XS Bundle

### DUAL FUNCTION OF WASHING AND CURING

- ❑ 360° all-round exposure
- ❑ Large washing volume of 7000ml
- ❑ Easy and intuitive operation
- ❑ UV Protection Cover
- ❑ Mercury XS Wash Machine Size: 255mm\*155mm\*385.5mm
- ❑ Maximum Wash Volume (with platform)180mm\*121mm\*153mm (without platform)201mm\*124mm\*255mm
- ❑ Mercury XS Cure Machine Size: 246mm\*230mm\*363.5mm



## Standard Photopolymer Beige & Grey Resin



### Features

- ❑ Low Shrinkage and High Precision
- ❑ Fast Curing and Great Stability
- ❑ Bright and Stunning Colours

### Specification

- ❑ Material: Polymer Resin
- ❑ Form: Liquid
- ❑ Compatible With: All LCD / MSLA / DLP 3D Printers
- ❑ Hardness: 84 D; Shrinkage: 7.1 %

## Dental Castable UV Resin

### Features

- ❑ High Precision
- ❑ Low contraction / viscosity
- ❑ Stable and fast healing

### Specification

- ❑ Solidification by wavelength: 385 to 410 nm
- ❑ Density: 1.05 - 1.25 g/cm<sup>3</sup>
- ❑ Shore hardness: 85 - 88 D
- ❑ Viscosity (25°C): 50 - 170 mPa.s



## Lab Brushless Motor

### OZ BLACK



#### DESKTOP TYPE STRONG TORQUE AT LOW SPEED VERTICAL / HORIZONTAL PLACEMENT

RPM	: Max 50,000 rpm
Torque	: 7.8 N. cm
Input	: 100-120, 220, 240V (50/60 Hz)
Output	: 230V

- Right & Left turning ability
- Hand/Foot speed control system
- Self-diagnosis, Error display, work memory

## Lab Brush Motor

### Strong 210/105L



#### NON STAGE SPEED SYSTEM WITH HIGH TORQUE UNIQUE DESIGN, HANDPIECE HOLDER

RPM	: Max 40,000 rpm
Torque	: 3 N.cm with 105L
Input	: 100, 110/320, 240V (50/60 Hz)
Output	: DC0-32V

- Right & Left turning ability
- Foot speed control system
- Foot speed pedal and foot on/off pedal options

## Lab Brushless Motor

# STRONG 204C/102L



### BASIC MODEL, MINI-COMPACT DESIGN LESS VIBRATION STANDARD CARBON BRUSH MOTOR


- RPM : Max 40,000 rpm
- Torque : 280 gf. cm
- Input : 100, 110/320, 240V  
(50/60 Hz)
- Output : DC0-32V
- Right & Left turning ability
- Foot on/off switch

## Surveyor



### MORE PRECISE WORK IS POSSIBLE WITH THE EXCLUSIVE HANDPIECE

- High Quality precision instrument for milling drilling setting attached models and tapping
- Smooth and micro-fit arm
- Surveyor holder using for Exclusive Handpiece (Strong 102LSII)




## LAB HAND PIECES




### Model: Forte 100 OIII

RPM : 50,000  
 Weight : 254 g  
 Dimensions : L-171, L1-81, L2-90, D-29, D1-19



### Model: Forte 100 IIP


RPM : 50,000  
 Weight : 225 g  
 Dimensions : L-165, L1-81, L2-84, D-28, D1-19



PRODUCT CODE : 02E7


### Model: Strong 105L

RPM : 40,000  
 Weight : 190 g  
 Dimensions : L-153, L1-80, L2-73, D-29, D1-15



### Model: Strong 102L

RPM : 40,000  
 Weight : 200 g  
 Dimensions : L-152, L1-79, L2-73, D-29, D1-15



## STAR ST WHITE ZIRCONIA

### STRENGTH AND TRANSLUCENCY

- STAR ST blanks are made from biocompatibility zirconiumdioxide.
- The yttrium oxide stabilizing projects the material against cracks and increase the tensile and compressive strength.
- The special grain size allocation inside the material and the added aluminum oxide also result in extra strength during milling and subsequently in clinical use.
- The outstanding mechanical characteristics, excellent chemical durability and the unbeaten biocompatibility.
- Combined with the translucent color predetermine STAR ST as the ideal material provider for dental milling systems.



**Translucency: 43%**

**Strength: 1250Mpa**



### Available Sizes

ST 98 x 10mm

ST 98 x 18mm

ST 98 x 12mm

ST 98 x 20mm

ST 98 x 14mm

ST 98 x 22mm

ST 98 x 16mm

ST 98 x 25mm

Items	ISO10993 Standards Requirements	Test results
Cytotoxicity test	Cytotoxicity <=1	No cytotoxicity
Sensitization test	Sensitization test not allowed	No sensitization
Oral mucosa irritation test	Irritation not allowed	No stimulation
Acute systemic toxicity test	Sensitization not allowed	No acute oral toxicity
Genotoxicity test	Avirulent	AMES test: negative

## STAR ST PRESHADED ZIRCONIA

### STAR ST HIGH TRANSLUCENCY, APPLY TO ANTERIOR & POSTERIOR TEETH, FULL CONTOUR AND ALL KINDS OF BRIDGES.

The ST High blank is non-toxic and non-radioactive with excellent biological properties, teeth color is like the real teeth with good permeability and aesthetics. It is with high strength, good toughness and excellent wear resistance, X-cera pre-shaded ST blanks is the raw material with excellent property for making crowns.



- Indications : Veneers, anterior crown/coping, bridge
- Shades : Shofu color Halo/Pro
- Sintering curve : Room temperature  $\dot{1}0^{\circ}\text{C}/\text{min}$  to  $900^{\circ}\text{C}$  (preserve for 30min)  $\dot{3}.5^{\circ}\text{C}/\text{min}$  to  $1530^{\circ}\text{C}$  (preserve for 2h)  $\dot{1}2^{\circ}\text{C}/\text{min}$  to  $800^{\circ}\text{C}$   $\dot{Z}$  follow the furnace cooling to room temperature

**Translucency: 43%**

**Strength: 1250Mpa**

### Available Sizes

ST 98 x 10mm

ST 98 x 18mm

ST 98 x 12mm

ST 98 x 20mm

ST 98 x 14mm

ST 98 x 22mm

ST 98 x 16mm

ST 98 x 25mm

### Ordering Information:

STAR ST PRESHADED ZIRCONIA

A1	A2	A3	A3.5
A4	B1	B2	C1
C2	C3	D2	D3

Items	ISO10993 Standards Requirements	Test results
Cytotoxicity test	Cytotoxicity $\leq 1$	No cytotoxicity
Sensitization test	Sensitization test not allowed	No sensitization
Oral mucosa irritation test	Irritation not allowed	No stimulation
Acute systemic toxicity test	Sensitization not allowed	No acute oral toxicity
Genotoxicity test	Avirulent	AMES test: negative

## STAR SHT MULTI LAYER ZIRCONIA BLANK

**SUPER-HIGH TRANSLUCENCY,  
APPLY TO ANTERIOR**

### AVAILABLE SIZES

SHT MULTI LAYER 98 x 10mm

SHT MULTI LAYER 98 x 12mm

SHT MULTI LAYER 98 x 14mm

SHT MULTI LAYER 98 x 16mm

SHT MULTI LAYER 98 x 18mm

SHT MULTI LAYER 98 x 20mm



### ORDERING INFORMATION: STAR SHT MULTI LAYER ZIRCONIA

A2

A3

Layer	Proportion	Vickers Hardness (HV)	Strength (MPa)	Color	Translucency
Layer-1	30%	≥1300	≥1050	Shallow	49%
Layer-2	15%			↓	
Layer-3	15%				
Layer-4	20%				
Layer-5	20%				



Projects	Detection Result
Sintering Density	$\leq 6.0\text{g/cm}^3$
Bending Strength	$\leq 1050\text{Mpa}$
Vickers Hardness	125HV
Chemical Solubility	$\leq 100\text{ug/m}^2$

## Regular providing pre-shaded blanks

Other shades in Vita Classical shades are available with MOQ request



## Indications



Full contour  
for anterior



Veneer



Full contour  
for posterior



Coping



Inlay/Onlay



Crown bridge 5-7 units



Crown bridge full arch



Post & core



Telescope



Maryland bridge

## Star Glass Ceramic

Star Glass Ceramic offers strong mechanical strength and durability, with high translucency, provide natural and lifelike esthetics.



Transmittance:	HT,LT,MT,MO
Transmittance:	Active concentration of uranium -238<,≤1.0Bq/g
Density after partial crystallization:	2.3-2.6g/cm <sup>3</sup>
Density after complete crystallization:	2.4-2.7g/cm <sup>3</sup>
Biaxial flexural strength:	≥300Mpa
Fracture toughness:	≥2Mpa.m <sup>1/2</sup>
Expansion coefficient:	10.5±0.5
Vickers hardness:	≥550HV2
Chemical solubility:	∅100ug/cm <sup>2</sup>

## AGKEM PMMA Dental Disc



PMMA Multilayer offers natural gradient esthetics, high color stability, and easy polishing.

### FEATURES

Product Name: PMMA Dental DISC

Brand: Agkem

Usage: CAD/CAM

Packaging: 1 Block per Box

Indications: Temporary crown & bridge restorations

### Available Sizes:-

- |                                  |                                  |
|----------------------------------|----------------------------------|
| <input type="checkbox"/> 98*10mm | <input type="checkbox"/> 98*18mm |
| <input type="checkbox"/> 98*12mm | <input type="checkbox"/> 98*22mm |
| <input type="checkbox"/> 98*14mm | <input type="checkbox"/> 98*25mm |
| <input type="checkbox"/> 98*16mm |                                  |

## AGKEM TITANIUM DISC TC 600

### FEATURES

Product Name: PMMA Dental DISC

Brand: Agkem

Material: Titanium

Surface Finish: Mill Finish

Usage: CAD/CAM


Process Type: Laser Cut

### Available Sizes:-

- |                                  |                                  |
|----------------------------------|----------------------------------|
| <input type="checkbox"/> 98*10mm | <input type="checkbox"/> 98*18mm |
| <input type="checkbox"/> 98*12mm | <input type="checkbox"/> 98*22mm |
| <input type="checkbox"/> 98*14mm | <input type="checkbox"/> 98*25mm |
| <input type="checkbox"/> 98*16mm |                                  |



## Milling Bur



### DC Bur For (ROLAND, VHF, ARUM, XT CERA) Milling Machine



**BUR HEAD SIZE**  
0.6, 1, 2 mm Head

### Also Available DLC BURS

Roland Milling Machine, VHF Milling Machine (K3, K4),  
Arum Milling Machine, XTCERA Milling Machine

### Features:

- ❑ Exceptionally sharp edge, preventing chipping & cracks of teeth
- ❑ Compatible with Zirconia
- ❑ Long life span
- ❑ Minimum risk of bur cracking due to the high tenacity
- ❑ Precision Milling Performance
- ❑ Balanced for Vibration-Free Operation
- ❑ Reduced Tool Change Frequency



## STAR POLVERE PREMIUM POWDER

Titanium powder dental alloy powder with improved properties for the production of dentures in the laser melting process

### Composition in percent by mass:

Ti	89%
Al	6%
V	4%
N,C,H,Fe,O	<1%

Packing: 5 kg

Size: - 10-30µm

## STAR POLVERE POWDER

Star polvere is a CoCr dental alloy powder with improved properties for the production of dentures in the laser melting process

### Composition in percent by mass:

Co	63.9%
Cr	24.7%
W	5.4%
Mo	5.0%
Si	1.0%

Packing: 5 kg

Size: - 10-45 µm



## STAR Ceram Micro Bond

- Ultra-thin protective layer (0.01-0.02 mm)
- Prevents edge darkening (garland area)
- Ceramic-based: ensures strong bond with opaquer
- Reusable casting cone
- Insulates metal oxides during firing
- Suitable for CoCr frames & low-noble alloys

Packing : 15 gm

## MESA MAGNUM CLARUM

- **One of the first produced by MESA**
- Approved by technicians around the world for over 35 years
- Ease of casting
- Excellent polishing and compactness
- Very high compatibility with various aesthetic coatings

### Composition in %:

Nickel (Ni) 63% Chrome (Cr) 25%

Molybdenum (Mo) 9%

Silicon (Si) 2%



## MESA-MAGNUM CERAMIC S

### One of the first produced by MESA

- Approved by technicians around the world for over 35 years
- Excellent mechanical properties
- Excellent polishing and compactness
- High compatibility with various aesthetic coatings

### Composition in %:

Nickel (Ni) 65% Chrome (Cr) 24%

Molybdenum (Mo) 10%

Silicon (Si) 1%





## STARBOND+

### without BERYLIUM

**NICKEL BASED DENTAL ALLOY FOR CERAMIC**

Composition	
Ni	64%
Cr	24%
Mo	10%
Si	1.5%
others	0.5%

## STARBOND+

### with BERYLIUM

**A HIGH STRENGTH NICKEL CHROME WITH BERYLIUM ALLOY FOR SUPERIOR BOND STRENGTH**

Composition	
Ni	76%
Cr	15%
Mo	5%
Be	1.7%
Al	2%
other	3% Max



# INVESTMENT MATERIALS

## SUNVEST

Dental Investment Material for Crown & Bridge

**Net Weight:** 400g / Pouch

**Carton Contains:** 30 Pouches & 3 ltr Liquid

**Mixture Ratio**

Powder : Liquid

100 gm : 24-25 ml



## STARVEST

Dental Investment Material for Partial Denture

**Net Weight:**

400g / Pouch

**Carton Contains:** 50 Pouches & 4 ltr Liquid

**Mixture Ratio**

Powder : Liquid

100 gm : 22-23 ml


## DENTROCK ULTRA HARD DIE STONE (TYPE-4)



- Very low expansion
- Splintering resistant
- Excellent casting under vibration
- Abrasion resistant
- High compressive strength

Dent Rock is an extra hard stone with high compressive strength, recommended for making master models for fixed prostheses.

Packaging:  
1 kg pack x 15 pcs.



## SNOW ROCK ULTRA HARD DIE STONE (TYPE-4)



### DIE STONE With Superior Quality

- Low setting expansion: 0.09%
- Glass-like accurate surface and beautiful color
- Excellent edge stability
- Available in color : Golden Brown, Green, Grey
- **Packaging**  
(1.5 kg Alu. pack x 12) - 18 kg

## DIE STONE EXTREME TYPE-4 (DIE STONE GOLD)



### SPECIAL FORMULA FOR IMPLANT & CAD/CAM

Super precision innovative Nano Technique silky surfaces with brilliant color ultra hard strength

- Low setting expansion: 0.07%
- Silky surfaces with brilliant color
- Good Scanability & Extremely Low roughness
- Innovative Nano Technique
- Available in color : Golden Brown
- **Packaging:** (1.5 kg Alu. pack x 8) - 12 kg




# CARBIDE BURS




High quality tungsten carbide lab burs with rapid cutting and smooth trimming of porcelain castings high performing burs with exceptionally long lifetimes. These include the unique soldering technique of carbide till head to shank unlike other burs in market with Carbide head and steel shank In addition, we achieve perfect concentricity resulting in smoother, vibration-free running.




Round Bur



Square Cylinder



Cross-Cut Square Cylinder



**Compatibility:** Smaller rotary tools.  
**Shank diameter:** 3mm(about 1/8").  
**Shank Length:** 35mm. From right to left direction of cutting uniform movement, not like the use of wheel for reciprocating cutting.

## CUTTING DISCS



### SEPARATING AND GRINDING DISK FOR METAL AND CERAMIC WORKING.



#### CUTTING DISCS

- ◆ Non-Pollutant
- ◆ Absolutely heatless
- ◆ For precise cutting of ceramics and precious alloy.
- ◆ Medium fine grit size

#### AVAILABLE SIZE:

Cutting Discs

(D32 x 2.2 mm),

(D38 x 0.7 mm),

(D38 x 1.0 mm)

Packing:100 pcs/box

## ARGON WHEELS & POINTS



Size

7/8 \* 1/4 inch



Size

7/8 \* 1/4 inch



Size

5/8 \* 1/16 inch



Size

7/8 \* 1/4 inch



Packaging: 100 pcs box  
Suggested Speed:30,000 rpm

◆ For polishing of semi precious alloys & metals of silicone carbide.

◆ For Pre-refing, Polishing, Smoothing

◆ Suitable for Chrome-Cobalt Alloy, semi-precious & titanium alloys

## Dental Procelain Furnace



- ◆ Double vacuum operation
- ◆ Crystallization
- ◆ Tempering over zirconia
- ◆ Auto vacuum calibration
- ◆ Max. Temperature : 1200°C
- ◆ Max. firing temperature : 1180°C
- ◆ Temperature ramp. : 1 to 120°C
- ◆ Programs : 200
- ◆ Acoustic signal : ✓
- ◆ Weight : 13.5 kg
- ◆ Width x Height x Depth (cm) : 23x57x24
- ◆ USB port connection : ✓
- ◆ Acoustic signal : ✓
- ◆ Cylindrical firing chamber : 95mm x 65mm

Solidilite V, the new generation indirect light curing unit features a compact, user friendly, sturdy yet ergonomic design with an integrated fan and a protective soft start cure made to optimize efficient light polymerization of indirect composite such as CERAMAGE.

## Solidilite 250V

- ◆ Dimensions : W 230mm x  
D 252mm x H 354mm
- ◆ Weight: approx : 9.5 kg
- ◆ Supply voltage : 200 V-240 V, 50/60 Hz
- ◆ Fuse : 2 x 3.15 A time-lag
- ◆ Lamp : 4 JCR 110V 150 W  
halogen lamps
- ◆ Spectral range : 400-550 nm
- ◆ Solidilite : 250 V



# VINTAGE PRO



## VINTAGE PRO BODY



### Ordering Information:

VINTAGE HALO BODY 50g

A1B	A2B	A3B	A3.5B	A4B	B1B
B2B	B3B	B4B	C1B	C2B	C3B
C4B	D2B	D3B	D4B		

VINTAGE HALO BODY 200g

A1B	A2B	A3B	A3.5B	A4B	B1B
B2B	B3B	B4B			

## VINTAGE PRO OPAQUE PASTE

### Ordering Information:

VINTAGE PRO Opaque Paste 6g

Base		A10	A20	A30	A3.50
A40	B10	B20	B30	B40	C10
C20	C30	C40	D20	D30	D40



## VINTAGE PRO BASE/SHADE OPAQUE POWDER



Ordering Information: VINTAGE PRO Paste opaque 50g

Base		A10	A20	A30	A3.50
A40	B10	B20	B30	B40	C10
C20	C30	C40	D20	D30	D40

## VINTAGE PRO INCISAL/STANDARD ENAMEL



A layer of translucent powder is recommended to be used along with standard enamel to give a more natural appearance.

### Ordering Information:

VINTAGE PRO Opal Enamel 50g and 200g

SHADE	57	58	59	60
-------	----	----	----	----



## VINTAGE PRO GUM



The PRO GUM porcelain offers an excellent range and armamentarium to best replicate gingival shades as close to the remaining natural gingiva thus enhancing the overall aesthetics of the restoration.

### Ordering Information:

VINTAGE PRO ADD-ON 50g

SHADE	LP	DP
-------	----	----

## VINTAGE MIXING LIQUID HC

The VINTAGE Mixing Liquid HC provides an ideal consistency and improves the overall handling properties of all porcelain powders.

### Ordering Information:

VINTAGE MIXING LIQUID HC

QTY	50 ML	500 ML
-----	-------	--------

**LIQUID SYSTEM**



## VINTAGE PRO POWDER OPAQUE LIQUID



- ◆ The VINTAGE PRO Powder Opaque Liquid is a water based medium which improves the overall workability when mixed along with the PRO opaque powder while masking the Inner metal Coping.

### Ordering Information:

VINTAGE PRO POWDER OPAQUE LIQUID

QTY	50 ML	500 ML
-----	-------	--------



## VINTAGE HALO BODY

# VINTAGE Halo



### Ordering Information:

VINTAGE HALO BODY 50g

A1B	A2B	A3B	A3.5B	A4B	B1B
B2B	B3B	B4B	C1B	C2B	C3B
C4B	D2B	D3B	D4B		

VINTAGE HALO BODY 200g

A1B	A2B	A3B	A3.5B	A4B	B1B
B2B	B3B	B4B			

## VINTAGE HALO INCISAL / STANDARD ENAMEL

### Ordering Information:

VINTAGE HALO INCISAL / STANDARD ENAMEL 50g & 200g

SHADE	57	58	59	60
-------	----	----	----	----



## VINTAGE HALO OPAQUE LIQUID



The VINTAGE Halo Opaque Liquid is a water based medium which improves the overall workability when mixed along with the Halo opaque powder while masking the inner metal coping.

### Ordering Information:

VINTAGE HALO OPAQUE LIQUID

QTY	50 ML	500 ML
-----	-------	--------

## VINTAGE HALO PAINT OPAQUE LIQUID



The VINTAGE Halo Paint Opaque Liquid is an oil based medium which creates an Ideal Paste like viscosity when mixed with Halo Opaque Powder thus enhancing the handling properties of the porcelain Powder.



### Ordering Information:

VINTAGE HALO PAINT OPAQUE LIQUID

QTY	50 ML	350 ML
-----	-------	--------

## VINTAGE AL/ZR OPAQUE LINER LIQUID



The AL/ZR opaque liner liquid is used to adjust the viscosity of the ZR Opaque Liner paste aiding uniform application over the coping in thin layers.

### Ordering Information:

VINTAGE AL/ZR OPAQUE LINER LIQUID

QTY	3 ML
-----	------

## VINTAGE MP BASE/SHADE OPAQUE PASTE

### Ordering Information:

VINTAGE MP OPAQUE PASTE 5g

Base		A10	A20	A30	A3.50
A40	B10	B20	B30	B40	C10
C20	C30	C40	D20	D30	D40





## VINTAGE ZR OPAQUE LINEAR PASTE

### Ordering Information:

VINTAGE ZR OPAQUE LINER PASTE 5g

A10	A20	A30	A3.50	A40	B10
B20	B30	B40	C10	C20	C30
C40	D20	D30	D40		



## VINTAGE ZR BODY / ENAMEL

### Ordering Information:

VINTAGE ZR BODY / ENAMEL 50g

A1B	A2B	A3B	A3.5B	A4B	B1B
B2B	B3B	B4B	C1B	C2B	C3B
C4B	D2B	D3B	D4B		

SHADE	OPAL 55	OPAL 56	OPAL 57
	OPAL 58	OPAL 59	OPAL 60



## VINTAGE ZR GUM

### Ordering Information:

VINTAGE ZR GUM

AVAILABLE IN 15g
GUM DP - DARK PINK
GUM LP - LITE PINK



## VINTAGE ART GLAZE PASTE



### Ordering Information:

VINTAGE ART  
GLAZE PASTE 3g

## VINTAGE ART SHADE STAINS

### Ordering Information:

VINTAGE ART SHADE STAINS

- ◆ AS - A Shade
- ◆ BS - B Shade
- ◆ CS - C Shade
- ◆ DS - D Shade



## VINTAGE ART BASE COLOUR STAINS



### Ordering Information:

VINTAGE ART BASE COLOR  
STAINS 3g

P - Pink	Y - Yellow	BL-Blue
O - Orange	V - Violet	G- Green
B - Black	W - White	G - Gray

## VINTAGE ART STAIN LIQUID

### Ordering Information:

VINTAGE ART STAIN LIQUID

**PACKAGING** 50 ML



## CERAMAGE POLISHING KIT



### Packaging :

The Kit Contains:-

- ◆ Robot Carbide fissure bur HP
- ◆ Silicone Points HR2 each of 1l. 13S & 2pc of 28
- ◆ SHOFU Pivot Brush x 1
- ◆ SHOFU Pivot Brush SC x 1
- ◆ SHOFU Mandrels x 2

## CERAMAGE MODELLING LIQUID

### Ordering Information:

VINTAGE MODELLING LIQUID

PACKAGING	6 ML	PIN	1991
-----------	------	-----	------



## CERAMAGE SPACER AND SEP



### Ordering Information:

CERAMAGE SPACER & SEP

PACKAGING	PIN
Spacer 7ml	1994
Separator 7ml	1993

## CERAMAGE OXY BARRIER

### Ordering Information:

CERAMAGE OXY-BARRIER

PACKAGING	10 ML	PIN	1992
-----------	-------	-----	------



## CERAMAGE FLOWABLE



Available in 12 Shades :

F-OD A3	F-A3B	F-ROOTAB	F-59
F-T-GLASS	F-CT	F-BT	F-BR
F-W	F-GUM-BR	F-GUM-V	F-GUM-R

Ordering Information: CERAMAGE FLOWABLE

## CERAMAGE GUM

Available in 3 types :

- GUM - L : Lite Pink
- GUM - D : Dark Pink
- GUM - OR : Orange



Ordering Information:  
CERAMAGE GUM

## VINTAGE SHADE GUIDE



Packaging:

Available as porcelain  
& resin type.

## PROCELAIN BRUSHES

Ordering Information:  
PORCELAIN BRUSHES

Packaging:

Porcelain Brush Set contains Brushes x 6 (No. OP0,  
No.P3, No.P5, No.P7, No.S0, No.S1)



## REFILL BRUSHES ARE INDIVIDUALLY PACKED



Packaging:

No. P-5	No. OP-0
No. P-7	No. S1
No. P3	



### Branch Offices:-

#### AHMEDABAD OFFICE

Shade No 74, Sankalp industrial state, Nr. S.P. Ring Road Nana Chiloda Ahmedabad-382330

#### BENGALURU OFFICE

Shop No:03, Ground Floor, Near Global IT Park, Pattanagere Main Road, Rajarajeshwari Nagara Bengaluru-560098, Karnataka

#### CHANDIGARH OFFICE

Shop No. 31, (Ground Floor) BhabatRoad Zirakpur, Sub-Tehsil Zirakpur, Punjab, - 140603

#### CHENNAI OFFICE

Plot no -41, Adinath Trade Complex, Ring road housing sector, Chennai, Tamil Nadu, 600060

#### GREATER NOIDA OFFICE

Site no - 4, d-11, kansa industrial area, greater Noida, Gautam Buddha Nagar, Uttar Pradesh, 201308

#### HYDERABAD OFFICE

Shop 6, Plot no. 44, Mamta Nagar Colony , Road No. 1, Nagole Hyderabad-68

#### INDORE OFFICE

Shop no-14, CM-2, Amar Aangan Basant Vihar, Indore, Madhya Pradesh - 452018

#### JAIPUR OFFICE

Shop no. 13-14 near Geetanjali Hotel Toll Tax, Tonk Road, Sanganer, Jaipur Ph. No.:0141 486 6081

#### KOLKATA OFFICE

Kolkata Office Address: 34A, Raja naba krishna street, Sovabazar, Kolkata-700005

#### MUMBAI OFFICE

412, 4th Floor Plot No.X21, MIDC Cross Road 21, Marol Industrial Area, Andheri East, Mumbai- 400069

## Agkem Impex Pvt. Ltd.

### Corporate Office

📍 503-504, Rishabh Corporate Tower, Plot no.16 Community Center, Karkardooma, Delhi-110092

☎ 011-43022244/55    ✉ info@agkem.com



Sales



Instagram



Youtube



Website



LinkedIn



Service/Support



Facebook